

Instructions for A.V. Fistula Needles

【 Scope of Application 】

- This product is suitable for establishing vascular access in hemodialysis treatment.

【 Main structure and properties 】

- The product is made up of protective cap, needle stand, tubing, on-off clamp, luer lock connector And luer lock endcap;
 - This product has a good property of rigidity, toughness and corrosion resistance;
 - Tubing is transparent that can easy to observe;
 - Tubing is flexible enough to resist winding and twisting;
 - On-off clamp can make the blood rate change from minimum to maximum;
 - The design of needle stand to make it convenient for handle.
-
- The diagram shows a medical device assembly. It includes a red protective cap, a grey needle, a clear tubing, a blue on-off clamp, a grey luer lock connector, and a grey luer lock endcap. The components are connected in a loop, with the tubing passing through the clamp and connecting the two luer lock ports. Labels with arrows point to each part: Protective cap, Needle, Luer lock endcap, Screw connector, Luer lock connector, On-off clamp, and Tubing.



【 Specification 】

- ・ 15G~17G

【 Usage 】

- Firstly, check expiry period and completeness of the products before using;
- Open the package from the peel site and take out A.V. fistula Needles;
- Connect connector with pipe of machine;
- Pull out needle protective cap and expel air in tubing;
- Push the needle spike into the body vein according aseptic manipulation and fix it.

【Precautions】

- It may causes hypothermia especially to weak patients;
- It may causes abnormal rhythm of ventricles of heart because of fast flow rate;
- There should be guardian during infusion in case of any accident;
- It should be sterile and free from pyrogen before package opened;
- It has a valid period of five years, Do not use if it is beyond valid date;
- Discard after single use.

【Warnings and Indications】

- Do not use if the package has been opened or damaged, protective cap has been removed, there are missing parts, and if foreign matter is present.
- It should be used by special trained person, and the skin should be correctly sterilized.
- The products can not be reused. The repetition use may cause cross infection, even serious illnesses such as hepatitis, HIV and other infectious disease.

【Storage】

- Product should be kept in a condition that relative humidity is less than 80%, it also should be free from corrosive gas, cool, dry with good ventilation and clean.

【Symbols】

Unique device identifier	Batch code	Catalogue number	Contains or presence of phthalates	Sterilized using ethylene oxide	
Date of manufacture	Manufacturer	Do not re-use	Caution	Do not use if package is damaged	
Medical device	Use-by date	Do not re-sterilize	Fragile, handle with care	Temperature limit	2862
Keep away from sunlight	Keep dry	Non-pyrogenic	This way up	Humidity limitation	
Single sterile barrier system		Sterile	REV.: B/0 DATE: 2025-05-30		