

Tkmd® Instructions for Pre-Filled Normal Saline Flush Syringe

【Indications for Use】

- The Pre-Filled Normal Saline Flush Syringe is intended to be used only for the flushing of indwelling vascular access devices.

【Product Description】

- The pre-filled normal saline flush syringe is intended for use in flushing compatible intravenous administration sets and indwelling intravenous access devices. It is made up of barrel, plunger, piston, nozzle cap and 0.9% sodium chloride injection.
- The pre-filled normal saline flush syringe is a three-piece, single use syringe with a 6% (luer) connector pre-filled with 0.9% sodium chloride injection, and sealed with a tip cap.

【Product Structure】

- It is made up of barrel, plunger, piston, nozzle cap and 0.9% sodium chloride injection.

【Product Specification】

- 3 ml, 5ml , 10ml, 3ml(05), 3ml(10), 5ml(10)

【Sterilization Method】

- Moist heat sterilization.

【Shelf life】

- 3 years.

【Usage】

- The clinicians and nurses should follow the steps below to use the product.
- Step 1: Tear the package at the cut part and take out the pre-filled normal saline flush syringe.
- Step 2: Push the plunger upward to release the resistance between the piston and the barrel.
Note: During this step do not unscrew the nozzle cap.
- Step 3: Rotate and unscrew the nozzle cap with sterile manipulation.
- Step4: Connect the product to an appropriated Luer connect or device.
- Step 5: The pre-filled normal saline flush syringe upwardly and expel all the air.
- Step 6: Connect the product to the connector, valve, or needleless system, and flush according to the relevant principles and the recommendations of the indwelling catheter manufacturer.
- Step7: Used pre-filled normal saline flush syringe should be disposed of in accordance with the requirements of hospitals and environmental protection departments. For single use only. Do not reuse.

【Contraindications】

- N/A.

【Cautions】

- Does not contain natural latex.
- Do not use if the package is opened or damaged;
- Do not use if the pre-filled normal saline flush syringe is damaged and leakage;
- Do not use if the nozzle cap is not installed correctly or is apart;
- Do not use if the solution is discolored, turbid, precipitated or any form of suspended particulate matter by visual inspection;
- Do not resterilize;
- Check the expiry date of the package, do not use if it is beyond the expiry date;
- For single use only. Do not reuse. Discard all the unused remaining parts;
- Do not contact the solution with incompatible medicines. Please review the compatibility literature.
- The product is not used for intramuscular injection or subcutaneous injection.
- People who are allergic to NaCl should not use the Pre-filled Normal Saline Flush Syringe

Tkmd® Instructions for Pre-Filled Normal Saline Flush Syringe

- Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established

【Storage】

- Product should be kept in a condition the temperature is 15-30℃, product cannot be frozen.





















【Risks of reusing single use devices】

- This product is a disposable product, Reprocessing, resterilization and/or reuse of this product is prohibited, reuse may cause infection or fail to reach the intended use.

【Residual risks】

- Failure to strictly enforce the instructions may result: Occlusion, Leakage, Phlebitis

【Graphs and symbols used for medical device labels】

| | | | |
|---|---|---|-------------------------------|
|  | Single use |  | Temperature limit |
|  | Don't use if the package is damaged. |  | Keep dry |
|  | Consult instructions for use |  | Date of manufacture |
|  | Manufacturer |  | Use-by date |
|  | Batch code |  | Do not resterilize |
|  | Non-pyrogenic |  | Unique Device Identifier |
|  | Medical device |  | This way up |
|  | Catalogue number |  | Keep away from sunlight |
|  | Fragile, handle with care |  | Stacking weight limit by mass |
|  | Sterile fluid path that has been sterilized using steam |  | Humidity limitation |

Edition: E/0 Date: 2025-03-28

CE 2862



Anhui Tiankang Medical Technology Co., Ltd.
No.228 Weiyi Road, Economic Development Zone, Tianchang City, Anhui, China.
Shanghai International Holding Corp. GmbH (Europe)
Eiffelstrasse 80 20537 Hamburg, GERMANY Tel: 0049-40-2513175