



受控

Instructions of Extension Tubes

【Scope of Application】

- The product is applied in clinical infusion in the operation, there is need to be connected to infusion containers, a variety of drug infusion at the same time, the extension of the infusion pipeline.

【Main structure and properties】

- Tube is transparent, with normal or corrected vision to find bubbles in line;
- Tube is flexible enough to resist winding and twisting;
- The design of closure-piercing device to make it convenient for handle and usage;
- All interfaces should be cases, cases should not fall off naturally, and easy to dismantle.
- After the clamps hold the canula, should have no liquid flow, and should not damage the catheter in continuous use.

【Specification】

- 50463H, 50465H, 32866H.

【Usage】

- Firstly, check valid period and completeness of package before using;
- Tear the package at the cut part and take out the extension tubes;
- Extension tubes joint connection on the infusion vessel joint;
- According to the aseptic operation, start infusion or infusion drugs.
- Infusion operation must have the professionals nursing to implement and guardian.

【Precautions】

- It may causes hypothermia-esp, weak patients;
- It may causes abnormal thymem of ventricles of heart because of fast flow rate;
- There should maintain the guardianship during the infusion in case of any accident;
- It should be sterile and free from pyrogen before package opened;
- It has a valid period of three years, Do not use if it is beyond valid date.
- Discard after single use.

【Warnings and indications】

- Do not use if the package is opened or damaged, or protective cap is aparted, any part is missed, And foreign matter is found;
- It should be used by special trained person.
- The products can not be reused. The repetition use may cause cross infection, even serious illnesses such as hepatitis, HIV and other infectious disease.

【Storage】

- Product should be kept in a condition that relative humidity is less than 80%, it also should be free from corrosive gas, cool, dry with good ventilation and clean.

ISO 8536-9-CL-P

Unique device identifier	Batch code	Catalogue number	Caution	Sterilized using ethylene oxide
Date of manufacture	Manufacturer	LATEX FREE	Do not re-use	Do not use if package is damaged
Medical device	Use-by date	Do not re-sterilize	Fragile, handle with care	Temperature limit
Keep away from sunlight	Keep dry	Non-pyrogenic	This way up	Humidity limitation
Single sterile barrier system	Sterile	Consult Instructions For Use	REV.: B/0 DATE: 2025-05-30	