



INSTRUCTIONS FOR DISPOSABLE SYRINGE WITH SAFETY NEEDLE

Indications for use

Intended for use by health care professionals for general purpose fluid aspiration/ injection.

Main structure and properties

- This product is made up of barrel, plunger, piston and safety needle.
- The barrel should be clean and transparent enough to easy observe.
- Barrel and piston matches well and it has a good property of sliding, can easy to use.
- Plunger is easy convenient for handle and usage.

Specifications

- Three-part syringes with luer lock or luer slip tip including 1ml to 60ml;
- Safety needle including 16G to 30G.

Usage

- Remove the Syringe and Safety Needle from the package. (see figure 1)
- Move the safety sheath away from the needle and toward the syringe barrel. (see figure 2)
- Remove the needle cover. (see figure 3 and 4)
- Draw up medication and perform procedure. (see figure 5 and 6)
- Position sheath in preparation for device activation. Using a one-handed technique, push the tab forward with your finger or thumb so that the sheath is less than 90 degrees from the needle. (see figure 7)
- Activate sheath. Position the sheath approximately 45 degrees to flat surface Press down with a motion until a distinct audible click is heard and Visually confirm that the needle is fully engaged under the lock. (see figure 8)
- Dispose of used needles and materials following the policies and procedures of your facility as well as federal and local regulations for Sharps disposal. (see figure 9)



Figure 1

Figure 2

Figure 3

Figure 4

Figure 5

Figure 6

Figure 7

Figure 8

Figure 9

Contraindications and Precautions

- Before injection, expel all the air in barrel from the needle upwardly.
- It should be sterile and non-pyrogenic before package opened.
- After use, this product should be discarded in the specified medical waste collection box, shall not be reused, and should be treated according to the relevant provisions of medical waste disposal.
- It has a valid period of five years, check valid period and completeness of package before using, do not use if it is beyond valid date.

Warnings

- Do not use if the package is opened or damaged, or protective cap is aparted, and part is missed, and foreign matter if found.
- It should be used by special trained person, and the skin should be correctly sterilized;
- The products can not be reused. The repetition use may cause cross infection, even serious illnesses such as hepatitis, HIV and other infectious disease.

Storage

Product should be kept in a condition that relative humidity is less than 80%, it also should be free from corrosive gas, cool, dry with good ventilation and clean.

UN Class device identifier	Batch code	Expiry date number	Warning	Material
1. Medical device	2. Expiry date	3. Expiry date	4. Expiry date	5. Expiry date
6. Expiry date	7. Expiry date	8. Expiry date	9. Expiry date	10. Expiry date
11. Expiry date	12. Expiry date	13. Expiry date	14. Expiry date	15. Expiry date
16. Expiry date	17. Expiry date	18. Expiry date	19. Expiry date	20. Expiry date
21. Expiry date	22. Expiry date	23. Expiry date	24. Expiry date	25. Expiry date
26. Expiry date	27. Expiry date	28. Expiry date	29. Expiry date	30. Expiry date
31. Expiry date	32. Expiry date	33. Expiry date	34. Expiry date	35. Expiry date
36. Expiry date	37. Expiry date	38. Expiry date	39. Expiry date	40. Expiry date
41. Expiry date	42. Expiry date	43. Expiry date	44. Expiry date	45. Expiry date
46. Expiry date	47. Expiry date	48. Expiry date	49. Expiry date	50. Expiry date
51. Expiry date	52. Expiry date	53. Expiry date	54. Expiry date	55. Expiry date
56. Expiry date	57. Expiry date	58. Expiry date	59. Expiry date	60. Expiry date
61. Expiry date	62. Expiry date	63. Expiry date	64. Expiry date	65. Expiry date
66. Expiry date	67. Expiry date	68. Expiry date	69. Expiry date	70. Expiry date
71. Expiry date	72. Expiry date	73. Expiry date	74. Expiry date	75. Expiry date
76. Expiry date	77. Expiry date	78. Expiry date	79. Expiry date	80. Expiry date
81. Expiry date	82. Expiry date	83. Expiry date	84. Expiry date	85. Expiry date
86. Expiry date	87. Expiry date	88. Expiry date	89. Expiry date	90. Expiry date
91. Expiry date	92. Expiry date	93. Expiry date	94. Expiry date	95. Expiry date
96. Expiry date	97. Expiry date	98. Expiry date	99. Expiry date	100. Expiry date

Anhui Tiankang Medical Technology Co., Ltd.
228, Wuyi Road, Economic Development
Zone, Tianshang City, Anhui Province, P.R. China.

Shanghai International Holding Corp. GmbH(Europe)
Platz 90 20027 Hamburg, Germany
Tel:0049-40-2513175