



**INSTRUCTIONS FOR
DISPOSABLE SYRINGE WITH SAFETY NEEDLE**

INDICATIONS FOR USE

The disposable syringe with safety needle is specially designed as a tool for vein and hypodermic injection solution to human body. It suits users of different ages.

Main structure and properties

- This product is made up of barrel, plunger, piston and needle with safety structure.
- This barrel should be clean and transparent enough to easy observe.
- Barrel and piston matches well and it has a good property of sliding, can easy to use.
- Plunger is easy convenient for handle and usage.

SPECIFICATION

- TKSWSN01S: 1ml syringe with fixed safety needle
- TKSWSN02S: 1ml luer slip syringe with detachable safety needle
- TKSWSN03S: 1ml luer lock syringe with detachable safety needle
- Needle size: 22G, 23G, 24G, 25G, 26G, 27G, 28G, 29G, 30G

USAGE

- First, Check valid period and completeness of package before using;
- Tear the package at the cut part and take out the syringe with needle;
- Push and pull plunger till it is well lubricated;
- Check that the needle is tightly connected to the syringe(for detachable needle only);
- Remove protective cap of needle, then do as sterile manipulation requirements;
- After injection, turn safety shield forward until you feel the positive stop and hear an audible click, so as to make needle fully covered.

Warnings & Precautions

- Expel all the air in barrel from the needle upwardly, in case gas boit come into being if air be injected into human body blood vessels or tumor come into being if gas be injected into subcutaneous;
- It should be sterile and free from pyrogen before package opened;
- It has a valid period of five years, do not use if it is beyond valid date;
- Discard after single use;
- Do not use if the package is opened or damaged, or protective cap is aparted or foreign matter is found;
- It should be used by special trained person, and the skin should be correctly sterilized;
- The products can not be reused. The repetition use may cause cross infection, even serious illnesses such as hepatitis, HIV and other infectious disease.
- After use, this product should be discarded in the specified medical waste collection box, shall not be reused, and should be treated according to the relevant provisions of medical waste disposal;

STORAGE

Product should be kept in a condition that relative humidity is less than 80%, it also should be free from corrosive gas, cool, dry with good ventilation and clean.

SYMBOLS

Unique device identifier	Batch code	Catalogue number	Caution	Sterilized using ethylene oxide
Date of manufacture	Manufacturer	Consult instructions For Use	Do not re-use	Do not use if package is damaged
Medical device	Use-by date	Do not re-sterilize	Fragile, handle with care	Temperature limit
Keep away from sunlight	Keep dry	Non-pyrogenic	This way up	Humidity limitation
Single sterile barrier system	Sterile			
				REV.: A/0 EFFECTIVE DATE: 2025.09.16

Anhui Tiankang Medical Technology Co., Ltd.
No. 228, Weiyl Road, Economic Development
Zone, Tianchang City, Anhui Province, P.R. China.

Shanghai International Holding Corp. GmbH(Europe)
Eiffestrasse 80 20537 Hamburg, Germany
Tel:0049-40-2513175